




## GENESIS

State of the Art Cartridge  
Based Inert Gas Purification  
Technology

**NO REGENERATION**

**Do you have a glovebox with dated purifier technology, or no purifier at all?**

- Upgrade with VAC's GENESIS cost effective independent O<sub>2</sub> and/or H<sub>2</sub>O purifiers.
- These purifiers will maintain 1 ppm O<sub>2</sub> and H<sub>2</sub>O with no regeneration required by you!
- Upgrade your pressure control system with VAC's GP-1 GENESIS pressure controller.
- Simple plug in components can be installed in any inert glovebox, regardless of manufacturer (contact VAC for box leak rate specifications).
- Cost effective oxygen and moisture analyzers also available.
- Purchase all three as a kit, or as separate components.
- Cartridges last 6 months before replacement. Ask your sales representative about our cartridge exchange program.



Main Offices and Manufacturing Facilities  
4652 W. Rosecrans Ave, Hawthorne, CA 90250  
(310) 644-0255

East Coast Office (978) 561-1763  
www.vac-atm.com info@vac-atm.com

Vacuum Atmospheres Company  
US Patents 9,878,283; 10,549,004  
MADE IN THE USA

NQA-1 Compliant

ISO 9001:2015 Certified